


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 03/19/16
Time Processed : 16:37:17
Injection Volume : 0.2
Sample Name : Secretin, human
Sample ID : C7713BC020-1
Theoretical MW : 3039.50
Observed MW : 3039.00

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH