


# Mass Spectrum



Sample Information  
Month-Day/Year Processed : 01/24/21  
Time Processed : 6:01:37 PM  
Injection Volume : 0.4  
Sample Name : HA Peptide  
Cat No. : RP111735  
Lot No. : E2264321A  
Theoretical MW : 1102.15  
Observed MW : 1101.8

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

Equipment : GK11010007  
Interface Bias : +3.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50%H2O/50%MeOH