


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 01/01/18
Time Processed : 13:45:20
Injection Volume : 0.8
Sample Name : Neuropeptide Y, human
Sample ID : U0626CL270-1
Theoretical MW : 4271.68
Observed MW : 4270.50

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H2O/50% MeOH