


Mass Spectrum



Sample Information
 Month-Day Processed : 11/24/20
 Time Processed : 8:31:52 PM
 Injection Volume : 0.3
 Sample Name : PACAP (1-38), human, ovine, rat
 Cat No. : RP10335
 Lot No. : E4917320K
 Theoretical MW : 4534.26
 Observed MW : 4533.6

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250
 Block Temp : 200

Equipment : GK11010007
 Interface Bias : +4.5kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2ml/min
 B.conc : 50%H₂O/50%MeOH