


Mass Spectrum



Sample Information

Month-Day Processed : 12/22/20
Time Processed : 10:27:06 PM
Injection Volume : 0.4
Sample Name : HA Peptide
Cat No : RP11735
Lot No : E5831820L
Theoretical MW : 1102.15
Observed MW : 1101.3

Interface : ESI

Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment

Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H₂O/50% MeOH

Equipment : GK11010007