


# Mass Spectrum



Sample Information  
Month-Day Processed : 02/05/21  
Time Processed : 5:19:08 PM  
Injection Volume : 0.5  
Sample Name : HA Peptide  
Cat No : RP11735  
Lot No : E3764021B  
Theoretical MW : 1102.15  
Observed MW : 1102.8

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

Equipment : GK11010007  
Interface Bias : +4.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50% H<sub>2</sub>O/50% MeOH