


Mass Spectrum



Sample Information
 Month-Day Processed : 02/08/21
 Time Processed : 9:54:16 PM
 Injection Volume : 0.3
 Sample Name : Atrial Natriuretic Peptide (ANP) (1-28), human, porcine
 Cat No. : RP11927
 Lot No. : E3608321B
 Theoretical MW : 3080.50
 Observed MW : 3079.6

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250
 Block Temp : 200

Equipment : GK11010007
 Interface Bias : +3.5 kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2 ml/min
 B.conc : 50%H2O/50%MeOH