


# Mass Spectrum



## Sample Information

User : Gary  
 Month-Day Processed : 08/20/18  
 Time Processed : 19:53:50  
 Injection Volume : 0.2  
 Sample Name : PACAP (1-38), human, ovine, rat  
 Cat No. : RPI0335  
 Lot No. : PE6741808  
 Theoretical MW : 4534.26  
 Observed MW : 4534.4

Interface : ESI  
 Nebulizing Gas Flow : 1.5L/min  
 CDL Temp : 250°C  
 Block Temp : 200°C

Interface Bias : +4.5kV  
 Drying Gas Flow : 5 L/min  
 T.Flow : 0.2ml/min  
 B.conc : 50%H2O/50%MeOH