


# Mass Spectrum



## Sample Information

Acquired by : Gary  
Month-Day Processed : 11/05/17  
Time Processed : 14:53:21  
Injection Volume : 0.4  
Sample Name :  $\beta$ -Amyloid (25-35)  
Cat. No. : RP10008  
Lot. No. : PE6871710  
Theoretical MW : 1060.27  
Observed MW : 1060.10

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250°C  
Block Temp : 200°C

Interface Bias : +4.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50% $\text{H}_2\text{O}$ /50% $\text{MeOH}$