


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 08/16/19
Time Processed : 10:43:38 PM
Injection Volume : 0.4
Sample Name : Calcitonin Gene Related Peptide (CGRP) (8-37), rat
Cat No. : RP11090
Lot No. : E5443219H
Theoretical MW : 3127.60
Observed MW : 3127.6

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200