


# Mass Spectrum



## Sample Information

Month-Day Processed : 06/23/21  
Time Processed : 2:56:15 PM  
Injection Volume : 0.4  
Sample Name : Calcitonin Gene Related Peptide (CGRP), rat  
SCat No : RP11095  
Lot No : E8659621F  
Theoretical MW : 3806.25  
Observed MW : 3806.4

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

Equipment : GK11010007  
Interface Bias : +4.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50% H<sub>2</sub>O/50% MeOH