


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 10/14/17
Time Processed : 20:47:21
Injection Volume : 0.2
Sample Name : hgp100(25-33)
Cat No: RP20344
Lot No: PE3351710
Theoretical MW : 1155.30
Observed MW : 1155.30

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias
Drying Gas Flow
T.Flow
B.conc

: +4.5 kV
: 5 L/min
: 0.2 ml/min
: 50% H₂O/50% MeOH