


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 03/02/17
Time Processed : 7:43:22
Injection Volume : 0.2
Sample Name : Calcitonin Gene Related Peptide (CGRP) (8-37), rat
Cat No. : RPI1090
Lot No. : PE0041702
Theoretical MW : 3127.60
Observed MW : 3127.65

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH