


Mass Spectrum



Sample Information

Acquired by: Gary
 Month-Day Processed: 02/06/17
 Time Processed: 14:25:49
 Injection Volume: 0.2
 Sample Name: β -Amyloid (1-42), human
 Cat. No.: RP10017
 Lot. No.: PE3561612
 Theoretical MW: 4514.10
 Observed MW: 4514.25

Interface

Interface: ESI
 Nebulizing Gas Flow: 1.5L/min
 CDL Temp: 250°C
 Block Temp: 200°C

Interface Bias

Interface Bias: +4.5 kV
 Drying Gas Flow: 5 L/min
 T.Flow: 0.2 ml/min
 B.conc: 50%H2O/50%MeOH