


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 12/30/16
Time Processed : 10:16:25
Injection Volume : 1
Sample Name : β -Amyloid (1-42), human
Cat No. : RP10017
Lot No. : PE 4061612
Theoretical MW : 4514.10
Observed MW : 4513.75

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H₂O/50%MeOH