


Mass Spectrum



Sample Information

Acquired by : Gary
 Month-Day Processed : 01/07/18
 Time Processed : 17:13:41
 Injection Volume : 0.2
 Sample Name : β -Amyloid (1-42), human
 Cat No.: RP10017
 Lot No.: PE3651803
 Theoretical MW : 4514.06
 Observed MW : 4514.00

Interface

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250°C
 Block Temp : 200°C

Interface Bias

Interface Bias : +4.5 kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2 ml/min
 B.conc : 50%H₂O/50%MeOH