


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 07/30/18
Time Processed : 17:30:14
Injection Volume : 0.2
Sample Name : CGRP (8-37), rat
Cat No. : RP11090
Lot No. : PE1741807
Theoretical MW : 3127.60
Observed MW : 3127.6

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250 °C
Block Temp : 200 °C

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH