


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 07/02/20
Time Processed : 05:43:20 PM
Injection Volume : 0.3
Sample Name : OVA G4 peptide
Cat No. : RP20398
Lot No. : E4912220F
Theoretical MW : 906.08
Observed MW : 905.6

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200