


# Mass Spectrum



## Sample Information

Month-Day Processed : 03/31/21  
 Time Processed : 10:41:21 PM  
 Injection Volume : 0.4  
 Sample Name : MDYKDHGDYKDHIDYKDDDDDK peptide  
 Cat No. : RP21087  
 Lot No. : E2681421C  
 Theoretical MW : 2861.87  
 Observed MW : 2861.2

Interface : ESI  
 Nebulizing Gas Flow : 1.5L/min  
 CDL Temp : 250  
 Block Temp : 200

Equipment : GK11010007  
 Interface Bias : +3.5 kV  
 Drying Gas Flow : 5 L/min  
 T.Flow : 0.2 ml/min  
 B.conc : 50%H<sub>2</sub>O/50%MeOH