


Mass Spectrum



Sample Information

Month-Day Processed : 05/08/21
Time Processed : 9:11:12 AM
Injection Volume : 0.4
Sample Name : Secretin, human
Cat No. : RP10249
Lot No. : E2664421E
Theoretical MW : 3039.50
Observed MW : 3038.8

Interface

:ESI
Nebulizing Gas Flow :1.5L/min
CDL Temp :250
Block Temp :200

Equipment

:GK11010007
Interface Bias :+3.5 kV
Drying Gas Flow :5 L/min
T.Flow :0.2 ml/min
B.conc :50%H₂O/50%MeOH