


Mass Spectrum



Sample Information		Interface	: ESI
Acquired by :	Gary	Interface Bias	: +4.5 kV
Month-Day Processed :	12/24/17	Drying Gas Flow	: 5 L/min
Time Processed :	12:48:34	T.Flow	: 0.2 ml/min
Injection Volume :	0.2	B. conc	: 50% H ₂ O/50% MeOH
Sample Name :	β-Amyloid (1-40)		
Cat. No. :	RP10004		
Lot. No. :	PE7921712		
Theoretical MW :	4329.82		
Observed MW :	4330.60		