


# Mass Spectrum



## Sample Information

User : Gary  
Month-Day Processed : 12/06/17  
Time Processed : 10:43:35  
Injection Volume : 0.2  
Sample Name : C-Peptide (3-33), human  
Cat No. : RPI1116  
Lot No. : PE2691712  
Theoretical MW : 3020.26  
Observed MW : 3021.15

## Interface

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250°C  
Block Temp : 200°C

## Interface Bias

Interface Bias : +4.5kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2ml/min  
B.conc : 50%H2O/50%MeOH