

Sample Name :  $\beta$ -Amyloid (25-35)  
 Cat. No. : RP10008  
 Lot. No. : PE1891707

Pump A : 0.065% trifluoroacetic in 100% water (v/v)  
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)  
 Total Flow: 1 ml/min  
 Wavelength: 220 nm


Time	Unit	Command	Value	Comment
0.01	Pumps	Pump A B.Conc	5	
25.00	Pumps	Pump A B.Conc	65	
25.01	Pumps	Pump A B.Conc	95	
31.00	Pumps	Pump A B.Conc	95	
31.01	Pumps	Pump A B.Conc	5	
40.00	Pumps	Pump A B.Conc	5	
40.01	Controller	Stop		

<<Column Performance>>

<Detector A>

Column : Inertsil ODS-3 4.6 x 250 mm

Chromatogram



1 Det.A Ch1 / 220nm

Peak Table

Detector A Ch1 220nm

Peak#	Ret. Time	Area	Height	Area %
1	13.888	3475	674	0.258
2	15.440	24700	1572	1.831
3	16.481	3118	531	0.231
4	17.177	766	169	0.057
5	17.531	1281717	196850	95.039
6	17.782	33327	5286	2.471
7	21.670	1522	251	0.113
Total		1348625	205334	100.000