


# Mass Spectrum



## Sample Information

Acquired by : Gary  
 Month-Day Processed : 07/25/18  
 Time Processed : 9:16:36  
 Injection Volume : 0.2  
 Sample Name : Calcitonin Gene Related Peptide (CGRP) (8-37), rat  
 Cat. No. : RP11090  
 Lot. No. : PE2211807  
 Theoretical MW : 3127.60  
 Observed MW : 3128.0

Interface : ESI  
 Nebulizing Gas Flow : 1.5L/min  
 CDL Temp : 250°C  
 Block Temp : 200°C

Interface Bias : +4.5 kV  
 Drying Gas Flow : 5 L/min  
 T.Flow : 0.2 ml/min  
 B.conc : 50% H2O/50% MeOH