


Mass Spectrum



Sample Information

Month-Day Processed : 12/03/20
 Time Processed : 11:59:52 PM
 Injection Volume : 0.3
 Sample Name : Calcitonin Gene Related Peptide (CGRP), rat
 Cat No. : RP11095
 Lot No. : E8644520L
 Theoretical MW : 3806.25
 Observed MW : 3805.2

Interface

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250
 Block Temp : 200

Equipment

Equipment : GK11010007
 Interface Bias : +4.5 kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2 ml/min
 B.conc : 50% H₂O/50% MeOH