


Mass Spectrum



Sample Information

Acquired by : Gary
 Month-Day Processed : 10/18/18
 Time Processed : 16:03:54
 Injection Volume : 0.2
 Sample Name : C-Peptide (3-33), human
 Cat. No. : RP11116
 Lot. No. : PE2941810
 Theoretical MW : 3020.28
 Observed MW : 3021.2

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250°C
 Block Temp : 200°C

Equipment : GK11010007
 Interface Bias : +4.5 kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2 ml/min
 B.conc : 50% H₂O/50% MeOH