


Mass Spectrum



Sample Information
Month-Day Processed : 08/18/21
Time Processed : 9:03:37
Injection Volume : 0.2
Sample Name : MDYKDHGDYKDHIDYKDDDDDK peptide
Cat No. : RP21087
Lot No. : E5493321H
Theoretical MW : 2861.87
Observed MW : 2861.6

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H₂O/50%MeOH