

Sample Name :C-Peptide (3-33), human
 Cat. No. : RP11116
 Lot. No. : PE7941902
 Time Processed :8:52:49
 Month-Day-Year Processed : 03/06/2019

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow:1 ml/min
 Wavelength:220 nm

Time	Unit	Command	Value	Comment
0.01	Pumps	Pump A B.Conc	5	
25.00	Pumps	Pump A B.Conc	65	
25.01	Pumps	Pump A B.Conc	95	
27.00	Pumps	Pump A B.Conc	95	
27.01	Pumps	Pump A B.Conc	5	
35.00	Pumps	Pump A B.Conc	5	
35.00	Controller	Stop		


<<Column Performance>>

<Detector A>

Column : Inertsil ODS-3 4.6 x 250 mm

Equipment: SS-CM-0286

Chromatogram



1 Det.A Ch1 / 220nm
 Detector A Ch1 220nm

Peak#	Ret. Time	Area	Height	Area %
1	16.780	1219	207	0.018
2	16.912	10513	1654	0.153
3	17.115	32116	5059	0.468
4	17.431	6577156	982116	95.771
5	17.813	234066	38496	3.408
6	21.679	9821	1399	0.143
7	24.020	2726	329	0.040
Total		6867617	1029262	100.000