


MS Spectrum Graph



Interface :ESI
 Nebulizing Gas Flow :1.5L/min
 CDL Temp :250°C
 Block Temp :200°C

Sample Information

Acquired by : Admin
 Date Acquired : 2017-1-13 14:38:44
 Injection Volume : 1
 Sample Name : c-Myc Peptide
 Cat No. : RP11731
 Lot No. : PE7211701
 Theoretical MW : 1203.30
 Observed MW : 1202.80

Interface Bias : +4.5 kV
 Drying Gas Flow :5 L/min
 T.Flow :0.2 ml/min
 B.conc :50%H2O/50%MeOH