


Mass Spectrum



Sample Information

Acquired by : Gary
 Month-Day Processed : 10/13/18
 Time Processed : 12:46:48
 Injection Volume : 0.2
 Sample Name : β -Amyloid (25-35)
 Cat. No : RP10008
 Lot. No. : PE1901810
 Theoretical MW : 1060.27
 Observed MW : 1060.2

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250°C
 Block Temp : 200°C

Equipment : GK11010007
 Interface Bias : +4.5 kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2 ml/min
 B.conc : 50% H₂O/50% MeOH